



MIRACLE involves 4 disciplines:



Planning & Navigation

- Planning in virtual reality
- Augmented reality in surgery
- Endoscope navigation



Medical Robotics & Mechatronics

- Robotic endoscope
- Telemanipulation and concurrent feedback
- Force sensing technologies



Smart Laser Osteotome

- Optimized laser technology
- Waveguides for high energy laser
- Real-time tissue monitoring and analysis



Smart Implants

- Personalized implants
- Smart implant design



WSS

WERNER SIEMENS - STIFTUNG

Project leaders:

Prof. Dr. Philippe Cattin philippe.cattin@unibas.ch

Prof. Dr. Dr. h.c. Hans-Florian Zeilhofer hf.zeilhofer@unibas.ch

Group leaders:

Prof. Dr. Georg Rauter georg.rauter@unibas.ch

Prof. Dr. Azhar Zam azhar.zam@unibas.ch

Project coordinator:

Dr. Constanze Pfeiffer constanze.pfeiffer@unibas.ch

MIRACLE Project

Department of Biomedical Engineering University of Basel Gewerbestrasse 14 4123 Allschwil Switzerland

Additional information:

www.miracle.dbe.unibas.ch





Department of Biomedical Engineering











Robot-Assisted Computer-guided LaserosteotomE (MIRACLE) combines natural sciences and medicine. It will enable minimally invasive laser osteotomies (bone cuts) using integrated miniaturized systems.



Planning & Navigation



Medical Robotics & Mechatronics



Smart Laser Osteotome



Smart Implants



The Planning and Navigation team develops novel navigation technology in order to control the robot-assisted laser system during surgery. One outcome is the groundbreaking SpectoVR software. SpectoVR is a renderer for three dimensional medical data such as Computer Tomography (CT) data sets and enables visualizing images in three dimensions using virtual reality. Since 2018 all operations with patients undergoing elective aneurysm repair at the University Hospital Basel, Switzerland are prepared with SpectoVR technology. In spring 2018, for the first time SpectoVR technology was used successfully during a consultation to inform a patient about his operation.

Group leader:


Prof. Dr. Philippe Cattin E-mail: philippe.cattin@unibas.ch



In order to achieve repeatable high precision cuts even in minimally invasive procedures, for the first time worldwide, the Bio-Inspired RObots for MEDicine-Lab (BIROMED-Lab) has developed a parallel-robotic endoscope-tip. This tip attaches to the target tissue, precisely guides the laser, and decouples mechanically from the endoscope and the robotic structure. Thus, if the patient is moved or if the surgeon touches the surgery robot, disturbances will not be transferred to the endoscope tip. Planned applications for this MIRACLE endoscope will be in the fields of orthopedics, craniomaxillofacial surgery, neurosurgery, otolaryngology, traumatology, and spinal column surgery.

Group leader:


Prof. Dr. Georg Rauter E-mail: georg.rauter@unibas.ch



The Biomedical Laser and Optics Group (BLOG) provides laser and optical solutions to build the laser for in-vivo laser osteotomy and deliver the laser beam inside the human body through an endoscope. The team develops a laser with feedback system which guarantees for extremely precise cuts of almost all shapes in minimally invasive surgery. The laser system will not only cut bone but also provide important information for the surgeons. It will be able to monitor the laser cutting process and to detect the type of tissue being cut in real-time. This enables the laser osteotome to immediately shut down in case of unexpected events like cutting through a nerve or other soft tissues that should be preserved.

Group leader:

Prof. Dr. Azhar Zam E-mail: azhar.zam@unibas.ch



The Smart Implants group focuses on developing a novel Minimal-Invasive Modular Implant System. The implant system consists of a plurality of customized implants that can be easily connected to each other, forming a preplanned and customized implant or bone plate arrangement. The design and shape of the customized implant or bone plate arrangement is individually defined by 3D software planning based on preoperative imaging. The connection between the customized implants can be done with many different connection possibilities, e.g. plates and screws connection or plug-in connection.

Group leader:

Prof. Dr. Dr. h.c. Hans-Florian Zeilhofer E-mail: hf.zeilhofer@unibas.ch

